

Biopsy procedure

Biopsies should be taken by a punch (diameter : 2mm) and placed in microtubes 1.5ml (model : MCT-150-C, Axygen) correctly labelled with the identification number of the mouse (Lab Id). Please exclusively use permanent black marker, do not use green or red marker.

Biopsies should be stored at -20°C before genotyping. They must be accompanied by a list including the name of the researcher, the date, the mouse line and the genotype of the parents (analyses request form).

The results are directly keyboarded in the main animal file.
For the mice that are not registered in the main file, please provide us an excel file (sent by e-mail to cig_genotyping@unil.ch) on which the results will be keyboarded.