


Background proteins in immunoprecipitations: what they are & how to reduce them



UNIL | Université de Lausanne

Patrice Waridel, Bastienne Jaccard, Alexandra Potts and Manfredo Quadroni
Protein Analysis Facility, Center for Integrative Genomics, University of Lausanne, Switzerland



LAUSANNE
PROTEIN ANALYSIS FACILITY

Introduction

Immunoprecipitation (IP) is a mainstream tool to study protein-protein interactions

Nonspecific background proteins occurring in IP's are a major problem, resulting in :

- Increased sample complexity, noise
- False positives
- False negatives

Goals

- 1) Characterize nonspecific background proteins in IP's
- 2) explore the origins of the background
- 3) find ways to reduce it

Conclusions :

- 1 - Survey of 10 mock IP's shows great variability and low predictive value but gives indications on the nature of the most recurrent background proteins
- 2 - A large portion of the background is due to protein binding to Sepharose beads
- 3 - Pelletable aggregates form during 4°C incubation of lysates. Pellet is larger if the lysate is treated with benzonase (RNase+DNase activity)
- 4 - Proteins found in the aggregate pellet contain most of the "frequent contaminants" found in IP's.
- 5 - A second centrifugation step before addition of beads drastically reduces level of background in IP's => **improved IP protocol !**

HYPOTHESIS : insoluble aggregates formed after lysis adhere to Sepharose beads and are responsible for a large part of the background. Formation of such aggregates is correlated with the presence of nucleic acids

Survey :

Identify and align protein lists from 10 negative control IP's

10 mock IP's { Varying : lab of origin, cell line, lysis method, IP protocols
Constant : human cells, Ab-based, use ProteinA-Sepharose
Constant : analytical pipeline by GeLC-MS , data analysis


Results :

1393 proteins identified combining all samples (1 pep. >95% , prot. >95% probability)

Great variability : only 15 proteins present in all samples

However : most ubiquitous proteins are most abundant ones (based on spectral counts)

Define arbitrarily : "frequent contaminants" detected in 7 or more samples : **115 prots**




Ribosomal proteins	30
Heterogeneous nuclear ribonucleoproteins	11
Keratins	7
Heat shock proteins	6
Elongation factors	6
Helicases	5
Nucleolar proteins	3
Replication factors	3
Other chaperones	3
Tubulins	2
Histones and associated	2
Exportins	2
Poly(RC)-binding proteins	2
ATP synthases	2
Actin and actin-bdg.	2
Chaperones	2
Others	29



Due to the large variability in the pool of recovered proteins, the data cannot be used to determine conclusively a list of contaminant proteins.

Frequent contaminants are constituted in majority by proteins from a few functional families.

Test the modified protocol (anti-PCNA IP) and MS analysis



ProtA-Sepharose	-	+	+	+
α-PCNA Ab	+	+	+	-
Cell lysate	-	+	+	+
Benzonase + spin	n/a	-	+	-



Most "frequent contaminants" are also found in pellet P2.


Aggregate removal in the **modified protocol (IP2)** leads to cleaner IP's

Which component generates the background ?


Sepharose CL-4B	+	-	-	-
ProteinA-Sepharose-CL4B	-	+	+	-
Human IgGs	-	-	+	+
Incubation w. cell lysate	+	+	+	-

Background appears to be due mainly to protein binding to the Sepharose CL-4B resin


* IgG light, heavy chain



Which step generates the background ?



IgGs	+	+	+	+	+
ProtA-Sepharose	+	+	-	-	-
Cell lysate	+	+	-	+	+
Benzonase	+	-	+	+	+
Spin before adding beads	+	-	-	+	+
Pellet before IP	-	-	-	+	+



Pelletable aggregates form during 4°C incubation of a 1% NP-40 cell lysate. A **modified protocol** was designed which included removal of these aggregates.

Materials & Methods :

IP's for the survey were from 6 different labs and performed in varying conditions. Common points were : 1) starting material was cultured human cells 2) IP was performed using antibodies against a protein or an epitope tag 3) Affinity resin was based on ProteinA-Sepharose.
IP's for method development: HeLa cells were lysed in 1% NP40, 1mM MgCl2, 150mM NaCl, protease inhibitors by pipetting and 5 min incubation at 4°C. Human IgGs were from Sigma-Aldrich. The anti-PCNA antibody was from Santa Cruz Biotechnology. 20 ul ProteinA-Sepharose CL-4B (GE Healthcare) slurry was used for immunoprecipitations. Benzonase (Merck) was used at a concentration of 125 U/ml lysate. Other steps were as described in the boxes.
MS Analysis : all proteins were eluted from the beads by boiling in 2%SDS sample buffer and subsequently separated along 2 cm on a 10% polyacrylamide gel. Whole gel lanes were excised from the gel into 5 fractions which were digested with trypsin. Dried extracts were analysed by nano-rp-LC-MS/MS on an LTQ-Orbitrap coupled to an Agilent nano 1100 HPLC system. Exported peptide MS/MS data were pooled for each lane and used for database searches in UNIPROT using Mascot 2.2 (10 ppm mass tolerance on precursors and 0.5 Da tolerance in MS/MS). Protein lists were validated and aligned with Scaffold 2.2 with the following thresholds : at least one peptide with min. 95% probability, minimum protein probability 95%